

# CLiC Microscopic Data Organization Database (Group 68)

Stone Chen, Yunwen Ji, Wenquan Li, Yuki Zhang

February 4, 2019

## 1 Requirement Analysis

### 1.1 Introduction

#### 1.1.1 Purpose

The database revolves around the CLiC raw data in building its connections with three most relevant components: the lab member, the experiment, and the analysis results. This architecture would benefit all users (most likely the lab members) in efficiently retrieving pertinent information, such as the individual responsible for a dataset, the experimental procedures, the file PATHs, etc.

#### 1.1.2 Scope and Special Design Choice

We limit the sphere of our interest to only the biomedical imaging component of the lab, as there is huge need for safe and efficient data management, particularly for this branch of the lab. Note that this database is not responsible for storing multidimensional arrays such as videos, images, and analysis related information; instead, those information are backed up in the lab cloud drive, of which the respective PATHs are saved as attributes to corresponding entities.

#### 1.1.3 Terminology

CLiC Microscopy: Convex Lense induced Confinement Microscopy

#### 1.1.4 Resources

<https://pubs.acs.org/doi/10.1021/ac101041s>. Please also see the appendix for a brief overview of this method.

## 1.2 Database Description

### 1.2.1 Entities and Attributes

- **Lab Members:** Those who are currently working or have worked in the lab, including research assistants and researchers. They both share the attribute name, contact information, and status. Name is key attribute as there is very low chance of duplication. Status can have value of active/inactive, indicating whether the individual still works here.
- **Research Assistant:** A research assistant (undergraduates) is a subclass of lab members, it has no additional attribute but unique relationships.
- **Researcher:** A researcher is a subclass of lab members with a special attribute named level, which indicates the degree/job level in the lab, for example, PhD. It also has unique relationship compared to research assistant
- **CLiC Raw Data:** The raw data is generated using a microscope installed with CLiC in the format of TIFF time series. This multidimensional array will not be saved directly in the database but rather backed up in the cloud drive, where the unique PATHs are captured by our database. This entity has an

artificial ID as key as well as an additional attribute date. Furthermore, since there are different types of data generated, which will subsequently be studied differently, this raw data contains two subclasses depending on the type of data.

- **Nanoparticle Data:** This subclass of raw data represent the data used in the studies of nanoparticle. It has one unique relationship compared to the other subclass.
- **DNA Unwinding Data:** This subclass of raw data represent the data used in the studies of DNA unwinding. It has one unique relationship compared to the other subclass.
- **Binding Analysis Result:** This entity represents the result analyzed by research assistant based on DNA unwinding data. These results are different from nanoparticle analysis results in the procedures performed as well as the information saved, which includes binding count, diffusion coefficient etc. Those results are not saved directly as attributes but backed up in the cloud drive, of which the PATH is stored as an attribute of this entity. In addition, this entity contains two other attribute: date and analysis type. Date is just the date those results are generated, while analysis type can have values of 'rough', 'decent' or 'finalized', which represents the degree of confidence in the accuracy of a particular result. We decide to save this as an attribute because this can help users decide whether or not they would like to get the PATH and check this result. For example, if a result is of analysis type 'rough', the user may not consider it valuable to its purposes. Finally, the key of this entity is an artificial ID.
- **Nanoparticle Analysis Result:** This entity represents the result analyzed by research assistant based on nanoparticle data. These results are different from DNA binding analysis results in the procedures performed as well as the information saved, which includes mean square displacement, number of particles etc. Those results are not saved directly as attributes but backed up in the cloud drive, of which the PATH is stored as an attribute of this entity. In addition, this entity contains two other attribute: date and MSD (Mean Square Displacement). Date is just the date those results are generated, while MSD can assume certain numerical values. We decided to save MSD as an attribute because again it is an important analysis result that a researcher would like to know before taking a closer look at the whole dataset. For instance, if a researcher decides to review all results with MSD larger than  $k$ , then this setup makes it convenient to do so. Finally, the key of this entity is an artificial ID.
- **CLiC Dashboard:** The entity represent the quality control result of CLiC raw data, which is performed before all analysis and presented in format of a dashboard. It includes general data information and specific quality information that includes the focus, SNR (signal-to-Noise Ratio) and contaminant level. All those information will be backed up in cloud drive, where the corresponding PATH is saved as an attribute of this entity. In addition, since quality information are the strict requirements to deem a raw data as valuable for analysis, we decided to save them as attributes so that users can quickly obtain these information before proceeding to analysis. Last but not least, this entity also has an attribute named date to store the date these results are generated. Finally, the key of this entity is an artificial ID.
- **Protocol:** This is the outline of experimental procedures, assays and safety information etc. Researchers formulate protocols before performing experiments so that they ensure they have made comprehensive planning. Each researcher has their own ways of naming their protocols, including typically their initials and some keywords regarding their experiment, which is always unique. Thus we use those protocol IDs as key of protocol. Since those protocols are usually saved and backed up somewhere, an attribute of their PATHs is also saved.
- **Microscope:** This entity represents the microscopes, which are used to perform experiments and generate raw data. There are currently two sets of microscopes in the lab and a third one is underway. Each device has a unique ID that starts with 'Nikon' and a number, for example 'Nikon 1'. Therefore, those unique microscope IDs are used as the key of this entity. Furthermore, those microscopes can be on temporary maintenance sometimes so an attribute named status is used; it assumes values of either 'active' or 'inactive'.

- **Experiment:** This weak entity represents an experiment that follows a certain protocol with a particular set of microscope. Often times, researchers replicate experiments using the same protocol to collect more data or verify their results. To represent this in ER model, experiment is designed to be a weak entity that is dependent on a protocol; it has a partial key named 'try number' that indicates the replicate number. For example, if try number is equal to 3, then this particular experiment is the third replicate of a certain protocol. Obviously this try number is not unique unless it is tailored to a specific protocol that has a unique protocol ID.

### 1.2.2 Relationships

- **is a (from lab members):** *Lab members* have two different subclasses: *researchers* and *research assistants*, who conduct distinctive tasks in the lab.
- **supervise:** A *researcher* supervises a *research assistant*. This is a many-to-many relationship because a researcher can supervise more than one research assistant and one research assistant can be supervised by more than one researcher. There is a participation constraint on the research assistant, which requires that every research assistant must be supervised by at least one researcher.
- **formulate:** A *researcher* formulates a *protocol*. This is a many-to-one relationship because a researcher can formulate many protocols but a protocol can only be formulated by one researcher. There is also a participation constraint on the protocol because every protocol must be formulated by at least one researcher, indicating that every protocol must be formulated by exactly one researcher incorporated with the key constraint.
- **generate:** A *researcher* generates a set of *CliC raw data* through an *experiment*. This is a many-one-many relationship because an experiment can only be conducted by one researcher but a researcher can conduct many experiments. In addition, a researcher can generate many sets of data and a set of data can be generated by multiple researchers.
- **using:** The *experiment* uses the *protocol*. This is a one-to-many relationship since one experiment can only apply one protocol while one protocol can be utilized by replicate experiments. There is a participation constraint on the experiment because every experiment must be conducted according to at least one protocol, indicating that every experiment uses exactly one protocol incorporated with the key constraint.
- **with:** An *experiment* is conducted with a *microscope*. This is a one-to-many relationship because an experiment can only use one microscope while one microscope can be used by plenty of experiments. There is also a participation constraint on the experiment because every experiment must be conducted with at least one microscope, indicating that every experiment must use exactly one microscope incorporated with the key constraint.
- **is a (CliC Raw Data):** The CliC new data has two subclasses: nanoparticle data and DNA unwinding data.
- **Assess:** The *research assistant* assesses the *CliC new data* to generate the *CliC dashboard*. This is a many-one-many relationship since one research assistant can assess more than one set of CliC raw data to get more than one CliC dashboard. One set of CliC raw data or one dashboard can only be assessed by one research assistant and one dashboard can only be generated from one set of raw data. There is a participation constraint on the CliC dashboard because every dashboard must be generated by at least one research assistant and at least one set of CliC raw data.
- **Analyze (DNA unwinding data or nanoparticle data):** The *research assistant* analyzes the *DNA unwinding data (nanoparticle data)* to generate *binding analysis results (nanoparticle analysis results)*. This is a many-one-many relationship since one research assistant can analyze more than one set of DNA unwinding data to generate more than one set of binding analysis results. One set of DNA unwinding data or binding analysis results can only be analyzed or obtained by one research assistant. One set of binding analysis results can only be produced by one set of DNA unwinding data. There is a participation

constraint on the binding analysis results since every set of binding analysis results must be achieved by at least one research assistant and at least one set of DNA unwinding data.

## 2 Functional Analysis

### 2.1 Overview

#### 2.1.1 Why do we need to store the CLiC data and its derivatives?

The lab members perform their research utilizing this microscopy technique extensively, which means that a huge amount of data is generated every day. Those data play a critical role in testing hypothesis and supporting conclusions in the day-to-day work of every researcher. Therefore, the lab has already set up cloud and local drives for backup purposes as well as a set of strict rules to follow in storing those data.

#### 2.1.2 Why do we need this database?

Despite those data are safely stored, researchers in the lab still complains about the difficulty to retrieve those data and their relevant information efficiently. For example, a researcher could not recall whether a dataset has been analyzed or not and spent a whole day looking for it on the cloud. He ended up finding nothing so he asked a research assistant to analyze it again which took him another whole day. This example can be even worse if that data has indeed been analyzed but just no where to find. Therefore, if there can be a magic button that one can simply press to find and present all they want in front of them, that will be exactly what we long to have. In practice, we can have something very similar which only requires more than one simple press — the database that store the links between CLiC data and all its relevant aspects!

### 2.2 Application Description

The central component of this database is the CLiC raw data and the application seeks to link all relevant components, which includes three major components: the lab members, the experiments, and the analysis results. Currently, those that need to be backed up on the hard drive are consisted of research protocol, CLiC raw data, CLiC dashboard, and CLiC data analysis results (in order of when they are generated). As soon as those are backed up on the cloud, their corresponding component in the database is expected to be updated accordingly. In this way, if one needs to retrieve relevant information on one dataset, all relationship that stems from the raw data entity will offer immediate feedback. Potentially, we would like the database to have the following functions:

#### 2.2.1 Experimental details? Is it legit?

Researchers and their PI usually review the experimental details before looking at the results to ensure the data is valuable to answer some questions. Given a particular dataset, this database can report on where the protocol is, which microscope was used as well as the number of replicates (try number) etc.

#### 2.2.2 Data analyzed? I need it NOW!!

The main intersection between the work of a researcher and a researcher assistant is one ultimate question — Is this dataset analyzed? This database will help them interact better and more efficiently since this question can be answered by checking whether there is a analysis result entity with date specified linked to the particular set of raw data. Researchers will be able to know whether their assistants are idling, while research assistants can have a clearer idea of the upcoming work to do.

#### 2.2.3 This dataset is of such poor quality!! Prove it.

From time to time, some dataset generated can be of poor quality and luckily we have quality control measures in place. This assessment can be performed by any research assistant and backed up, however more importantly, those who will potentially be assigned this analysis work need to be aware of it so that they do not waste their

time. Rest assure, our database can capture this information through the relationship between raw data and CLiC Dashboard. Some of the hardline quality control parameters are saved directly as attributes, so that in most cases, with a quick reference of those values, a dataset can be trashed immediately if not qualified.

#### 2.2.4 Data analyzed!! Great, where is the one I want?

The greatest pain often comes from actually finding the particular subsets of data out of the pools of dataset. For example, find PATHs of all DNA unwinding analysis result in after Sept, 2018 that are decent counts. This is such a tedious task that every researcher suffers from in the lab. Nevertheless, with this database, all the PATHs can be output immediately with a line of query and potentially we can implement this as a user-friendly function (the magic button).

#### 2.2.5 Who is responsible??

This is a general question asked for everything in the lab. Conveniently, our database captures the person responsible for every task, or in other way around, what a person is responsible for in his work. This even includes the supervision responsibility between researcher and research assistant.

### 3 Relations

#### Entities

labMember (Name, contactInformation, Status)

Researcher (Name, contactInformation, Status, Level)

researchAssistant (Name, contactInformation, Status)

CLiCDashboard (dID, Date, PATH, Focus, SNR, Contaminants, crID, Name) (crID ref CLiC Raw Data(crID), Name ref researchAssistant (Name))

nanoparticleAnalysisResults (narID, Date, PATH, MSD, Name, crID) (crID ref nanoparticleData (crID), Name ref researchAssistant (Name))

bindingAnalysisResults (barID, Date, PATH, nnalysisType, Name, crID) (crID ref DNAUnwindingData (crID), Name ref researchAssistant (Name))

CLiCRawData (crID, Date, PATH)

nanoparticleData (crID)

DNAUnwindingData (crID)

Protocol (ProtocalID, PATH, Name) (Name ref Researcher (Name))

Microscope (MicroscopeID, Status)


#### Weak Entity

Experiment(protocolID, tryNumber, microscopeID, crdID, Name) (microscopeID ref Microscope(microscopeID), protocolID ref Protocol(protocolID), crdID ref CliCRawData(crdID), Name ref Researcher(Name))

#### Relationship

Supervise (Period, rName, raName) (rName ref Researcher (Name), raName ref researchAssistant (Name))

We did not feel there was a way to combine any relations. For example, if we did not use a relation table for the ER relationship **Supervise**, but save this as an attribute in the **researchAssistant** table, we would not be able store past records. This sitation could occur if a research assistant had worked with multiple researcher.



# CLiC Microscopic Data Organization Database P2 (Group 68)

Stone Chen, Yunwen Ji, Wenquan Li, Yuki Zhang

March 2, 2019

## 1 Modified Relations

The feedback we received is as follows:

*Sub-entities in ISA hierarchies should explicitly reference the super-entity in the relational translation.*

### Entities

LabMember (Name, contactInformation, Status)

Researcher (Name, contactInformation, Status, Level) (Name ref LabMember(Name))

ResearchAssistant (Name, contactInformation, Status) (Name ref LabMember(Name))

CLiCDashboard (dID, Date, PATH, Focus, SNR, Contaminants, crID, Name) (crID ref CLiC Raw Data(crID), Name ref researchAssistant (Name))

nanoparticleAnalysisResults (narID, Date, PATH, AvgRadius, Name, crID) (crID ref nanoparticleData (crID), Name ref researchAssistant (Name))

bindingAnalysisResults (barID, Date, PATH, analysisType, Name, crID) (crID ref DNAUnwindingData (crID), Name ref researchAssistant (Name))

CLiCRawData (crID, Date, PATH)

nanoparticleData (crID)(crID ref CLiCRawData(crID))

DNAUnwindingData (crID)(crID ref CLiCRawData(crID))

Protocol (ProtocolID, PATH, Name) (Name ref Researcher (Name))

Microscope (MicroscopeID, Status)

### Weak Entity


Experiment(protocolID, tryNumber, microscopeID, crdID, Name) (microscopeID ref Microscope(microscopeID), protocolID ref Protocol(protocolID), crdID ref CliCRawData(crdID), Name ref Researcher(Name))

### Relationship

Supervise (rName, raName, StartDate, EndDate) (rName ref Researcher (Name), raName ref researchAssistant (Name))

## 2 Database Construction

Please refer to ...sql for more information. The tables created after running all the setup commands are demonstrated in Fig.1.



Powered by yFiles

Figure 1: Database framework after running all create table commands. This representation is generated with DBVisualizer. Bold Text = "NotNull"; Key sign = "Primary Key"; Blue Arrow = "Foreign Key"

### 3 Manual Insert

Please refer to *q3.sql* for details on the insert commands. We inserted multiple triplet data into the *labmember* relation as well as its sub-classes *researchassistant* and *researcher*. In particular, *name* and *status* contain real data, while *contactinformation* is randomly generated. The newly inserted information in each of the three tables is shown, respectively:

```
[db2 => SELECT * FROM labmember
[db2 (cont.) => ;

          NAME           CONTACTINFORMATION STATUS
-----+-----+-----+-----+
Wendy Ji          (514)-561-5587    active
Stone Chen        (514)-531-3322    active
Zhi Zhang         (514)-557-3232    active
Zach Friedenberger (514)-324-1237    active
Gracia Gu         (514)-545-5123    active
Zhi Yue Zhang     (514)-512-5567    active
Cindy Shaheen     (514)-532-8521    active
Francis Stabile   (514)-274-4887    active
Daniel Berard      (514)-244-9817    active
Kim Metera         (514)-274-4921    active
Radin Tahvildari (514)-284-1237    active
Romain Berti       (514)-224-9387    active
Sabrina Leslie     (514)-224-2837    active
Shane Scott         (514)-274-9287    inactive
Haoran Liao         (514)-271-1927    inactive
Yash Patel          (514)-292-1927    inactive

16 record(s) selected.

[db2 => SELECT * FROM researcher;

          NAME           LEVELS
-----+-----+
Zhi Yue Zhang      master
Cindy Shaheen      PhD
Francis Stabile    PhD
Daniel Berard       associate
Kim Metera         associate
Radin Tahvildari  associate
Romain Berti       associate
Sabrina Leslie      PI
Shane Scott          PhD

9 record(s) selected.

[db2 => SELECT * FROM researchassistant
[db2 (cont.) => ;

          NAME
-----+
Gracia Gu
Haoran Liao
Stone Chen
Wendy Ji
Yash Patel
Zach Friedenberger
Zhi Zhang

7 record(s) selected.
```

## 4 Auto Data Generation

Please refer to *autoDataGen.py* and *autoDataGen.sql* for details on the automatically generated *Insert* commands. We used python to automatically generate 200 *Insert* commands for each table. Their connections are randomly constructed but still stick to the rules specified for each relation. Some inserted information are shown below:

```
db2 => SELECT * FROM protocol FETCH FIRST 10 ROWS ONLY
PROTOCOLID PATHS
-----NAME
2 C:\Protocol\2.pdf Kim Metera
3 C:\Protocol\3.pdf Cindy Shaheen
4 C:\Protocol\4.pdf Radin Tahvildari
5 C:\Protocol\5.pdf Radin Tahvildari
6 C:\Protocol\6.pdf Francis Stabile
7 C:\Protocol\7.pdf Radin Tahvildari
8 C:\Protocol\8.pdf Daniel Berard
9 C:\Protocol\9.pdf Francis Stabile
10 C:\Protocol\10.pdf Cindy Shaheen
11 C:\Protocol\11.pdf Daniel Berard

10 record(s) selected.

db2 => SELECT * FROM experiment FETCH FIRST 10 ROWS ONLY
PROTOCOLID TRYNUMBER MICROSCOPEID CRID NAME
----- -----
2 1 Nikon01 183 Daniel Berard
2 2 Nikon01 11 Daniel Berard
2 3 Nikon01 33 Daniel Berard
2 4 Nikon01 50 Daniel Berard
3 1 Nikon01 33 Cindy Shaheen
3 2 Nikon02 175 Cindy Shaheen
4 1 Nikon02 35 Sabrina Leslie
4 2 Nikon02 93 Sabrina Leslie
4 3 Nikon02 185 Sabrina Leslie
4 4 Nikon02 174 Sabrina Leslie

10 record(s) selected.

db2 => SELECT * FROM BINDINGANALYSISRESULTS FETCH FIRST 10 ROWS ONLY
BARID DATES PATHS ANALYSISTYPE CRID NAME
----- -----
1 12/16/2018 C:\BindingAnalysisResults\1 rough 154 Stone Chen
2 12/22/2018 C:\BindingAnalysisResults\2 rough 57 Zach Friedenberger
3 1/1/2019 C:\BindingAnalysisResults\3 rough 24 Stone Chen
4 1/8/2019 C:\BindingAnalysisResults\4 rough 167 Wendy Ji
5 6/01/2018 C:\BindingAnalysisResults\5 accurate 128 Stone Chen
6 6/12/2019 C:\BindingAnalysisResults\6 rough 95 Wendy Ji
7 8/05/2018 C:\BindingAnalysisResults\7 decent 188 Stone Chen
8 8/06/2018 C:\BindingAnalysisResults\8 accurate 156 Stone Chen
9 9/23/2018 C:\BindingAnalysisResults\9 rough 124 Zhi Zhang
10 12/18/2018 C:\BindingAnalysisResults\10 accurate 178 Zhi Zhang

10 record(s) selected.

db2 => SELECT * FROM CLICDASHBOARD FETCH FIRST 10 ROWS ONLY
DID DATES PATHS FOCUS SNR CONTAMINANTS CRID NAME
----- -----
1 06/27/2018 C:\CLICDashBoard\1 +1.50000000000000E+000 -5.47000000000000E+000 2539 175 Zhi Zhang
2 03/23/2019 C:\CLICDashBoard\2 +1.76000000000000E+000 +1.21000000000000E+000 600 186 Stone Chen
3 05/07/2018 C:\CLICDashBoard\3 +4.82000000000000E+000 +9.21000000000000E+000 145 43 Stone Chen
4 10/04/2018 C:\CLICDashBoard\4 +3.06000000000000E+000 +6.63000000000000E+000 2646 127 Zach Friedenberger
5 08/22/2018 C:\CLICDashBoard\5 +2.90000000000000E+000 +4.62000000000000E+000 1957 128 Zhi Zhang
6 07/03/2018 C:\CLICDashBoard\6 +2.95000000000000E+000 +9.03000000000000E+000 636 88 Stone Chen
7 07/04/2018 C:\CLICDashBoard\7 +4.77000000000000E+000 +1.31000000000000E+000 1845 192 Stone Chen
8 10/04/2018 C:\CLICDashBoard\8 +1.12000000000000E+000 +8.79000000000000E+000 3427 153 Zach Friedenberger
9 07/24/2018 C:\CLICDashBoard\9 +4.01000000000000E+000 +1.40000000000000E+000 773 1 Wendy Ji
10 08/28/2018 C:\CLICDashBoard\10 +3.09000000000000E+000 +5.41000000000000E+000 1288 115 Stone Chen

10 record(s) selected.

db2 => SELECT * FROM CLICDASHBOARD FETCH FIRST 10 ROWS ONLY
DID DATES PATHS FOCUS SNR CONTAMINANTS CRID NAME
----- -----
1 06/27/2018 C:\CLICDashBoard\1 +1.50000000000000E+000 -5.47000000000000E+000 2539 175 Zhi Zhang
2 02/23/2019 C:\CLICDashBoard\2 +1.76000000000000E+000 +1.21000000000000E+000 600 186 Stone Chen
3 03/23/2019 C:\CLICDashBoard\3 +4.82000000000000E+000 +9.21000000000000E+000 145 43 Stone Chen
4 10/04/2018 C:\CLICDashBoard\4 +3.06000000000000E+000 +6.63000000000000E+000 2646 127 Zach Friedenberger
5 08/22/2018 C:\CLICDashBoard\5 +2.90000000000000E+000 +4.62000000000000E+000 1957 128 Zhi Zhang
6 07/03/2018 C:\CLICDashBoard\6 +2.95000000000000E+000 +9.03000000000000E+000 636 88 Stone Chen
7 09/05/2018 C:\CLICDashBoard\7 +4.71000000000000E+000 +1.85000000000000E+000 1845 192 Stone Chen
8 10/06/2018 C:\CLICDashBoard\8 +1.12000000000000E+000 +8.89000000000000E+000 3427 153 Zach Friedenberger
9 07/24/2018 C:\CLICDashBoard\9 +4.01000000000000E+000 +1.40000000000000E+000 773 1 Wendy Ji
10 08/28/2018 C:\CLICDashBoard\10 +3.09000000000000E+000 +5.41000000000000E+000 1288 115 Stone Chen

10 record(s) selected.

db2 => SELECT * FROM CLICRANCDATA FETCH FIRST 10 ROWS ONLY
CRID DATES PATHS
----- -----
1 07/29/2018 C:\CLICRawData\1.mat
2 11/15/2018 C:\CLICRawData\2.mat
3 07/17/2018 C:\CLICRawData\3.mat
4 07/11/2018 C:\CLICRawData\4.mat
5 05/07/2018 C:\CLICRawData\5.mat
6 10/20/2018 C:\CLICRawData\6.mat
7 10/20/2018 C:\CLICRawData\7.mat
8 05/14/2018 C:\CLICRawData\8.mat
9 06/11/2018 C:\CLICRawData\9.mat
10 07/13/2018 C:\CLICRawData\10.mat

10 record(s) selected.

db2 => SELECT * FROM MICROSCOPE FETCH FIRST 10 ROWS ONLY
MICROSCOPEID STATUS
----- -----
Nikon01 active
Nikon02 active
Nikon03 inactive

3 record(s) selected.

db2 => SELECT * FROM NANOPARTICLEANALYSISRESULTS FETCH FIRST 10 ROWS ONLY
NARID DATES PATHS AVERADIUS CRID NAME
----- -----
0 02/18/2019 C:\NanoparticleAnalysisResults\1 +4.88000000000000E+001 85 Stone Chen
1 12/26/2018 C:\NanoparticleAnalysisResults\2 +3.70000000000000E+001 6 Wendy Ji
2 08/08/2018 C:\NanoparticleAnalysisResults\3 +4.17500000000000E+001 102 Wendy Ji
3 08/20/2018 C:\NanoparticleAnalysisResults\4 +4.21600000000000E+001 193 Stone Chen
4 08/09/2018 C:\NanoparticleAnalysisResults\5 +3.45600000000000E+001 152 Zhi Zhang
5 08/09/2018 C:\NanoparticleAnalysisResults\6 +3.42000000000000E+001 189 Zach Friedenberger
6 08/07/2018 C:\NanoparticleAnalysisResults\7 +3.64000000000000E+001 189 Zach Friedenberger
7 08/06/2018 C:\NanoparticleAnalysisResults\8 +2.87700000000000E+001 29 Stone Chen
8 05/20/2018 C:\NanoparticleAnalysisResults\9 +2.25200000000000E+001 64 Stone Chen
9 12/28/2018 C:\NanoparticleAnalysisResults\10 +4.97600000000000E+001 133 Zach Friedenberger

10 record(s) selected.
```

```

db2 => SELECT * FROM NANOPARTICLEDATA FETCH FIRST 10 ROWS ONLY
CRID
-----
1
2
3
4
5
6
7
8
9
10

10 record(s) selected.

db2 => SELECT * FROM SUPERVISE FETCH FIRST 10 ROWS ONLY
STARTDATE ENDDATE RNAME RNAME
----- -----
05/01/2018 - Radin Tahvildari Wendy Ji
05/01/2018 - Radin Tahvildari Zach Friedenberger
02/01/2017 06/01/2018 Radin Tahvildari Haoran Liao
02/01/2017 06/01/2018 Radin Tahvildari Yash Patel
05/01/2018 - Cindy Shaheen Stone Chen
05/01/2018 - Cindy Shaheen Zhi Zhang
05/01/2018 - Cindy Shaheen Gracia Gu

7 record(s) selected.

db2 => SELECT * FROM RESEARCHASSISTANT FETCH FIRST 10 ROWS ONLY
NAME
-----
Gracia Gu
Haoran Liao
Stone Chen
Wendy Ji
Yash Patel
Zach Friedenberger
Zhi Zhang

7 record(s) selected.

db2 => SELECT * FROM RESEARCHER FETCH FIRST 10 ROWS ONLY
NAME LEVELS
----- -----
Zhi Yue Zhang master
Cindy Shaheen PhD
Francis Stabile PhD
Daniel Berard associate
Kim Metera associate
Radin Tahvildari associate
Romain Berti associate
Sabrina Leslie PI
Shane Scott PhD

9 record(s) selected.

```

## 5 Analytical Query

### 5.1 Digging around!!

List the paths of all data analysis results (including binding analysis and nano-particle analysis) analyzed by an active researcher named "Wendy Ji":

```

1 SELECT PATHS
2 FROM NANOPARTICLEANALYSISRESULTS
3 WHERE NAME = 'Wendy Ji'
4 UNION
5 SELECT PATHS
6 FROM BINDINGANALYSISRESULTS
7 WHERE NAME = 'Wendy Ji';
7:25 (144) INS
Log 1: NANOPARTICLEANALYSISRESULTS (93) x
* PATHS
22 C:\BindingAnalysisResults\179
23 C:\BindingAnalysisResults\18
24 C:\BindingAnalysisResults\182
25 C:\BindingAnalysisResults\185
26 C:\BindingAnalysisResults\186
27 C:\BindingAnalysisResults\187
28 C:\BindingAnalysisResults\18
29 C:\BindingAnalysisResults\199
30 C:\BindingAnalysisResults\23
31 C:\BindingAnalysisResults\25
32 C:\BindingAnalysisResults\27
33 C:\BindingAnalysisResults\28
34 C:\BindingAnalysisResults\4
35 C:\BindingAnalysisResults\44
36 C:\BindingAnalysisResults\47
37 C:\BindingAnalysisResults\50
38 C:\BindingAnalysisResults\6
39 C:\BindingAnalysisResults\77
40 C:\BindingAnalysisResults\93
41 C:\BindingAnalysisResults\94
42 C:\BindingAnalysisResults\95
43 C:\NanoparticleAnalysisResults\101
44 C:\NanoparticleAnalysisResults\105
45 C:\NanoparticleAnalysisResults\114
46 C:\NanoparticleAnalysisResults\118
Format: Text
Time Status Command Exec Fetch Rows Message SQL/Command
▼ 22:25:53 STARTED Executing for: 'cs421g' [DB2 LUW], Schema: CS421G68
→ 22:25:53 SUCCESS SELECT 0.025 0.018 93 Result set fetched
▲ 22:25:54 FINISHED 0.025 0.018 93 Success: 1

```

### 5.2 Who has been slacking?

Rank active research assistants (name) by the number of CLiC raw data they have analyzed during January and February of 2019:

```

9
10 SELECT C.NAME, COUNT(UNIQUE(CRID)) AS NUMBEROFRDATA
11 FROM CLICDASHBOARD C, LABMEMBER L
12 WHERE C.DATES >= '2019-01-01' AND C.DATES <= '2019-02-28' AND C.NAME = L.NAME AND L.STATUS='active'
13 GROUP BY C.NAME
14 ORDER BY NUMBEROFRDATA;
15

```

7:25 [144] [INS]

Log 1: CLICDASHBOARD [4] X

SQL/Command

NAME	NUMBEROFRDATA
Zhi Zhang	4
Zach Friedenberger	8
Stone Chen	12
Wendy Ji	12

Time Status Command Exec Fetch Rows Message SQL/Command

22:28:51 STARTED Executing for: 'cs421g' [DB2 LUW], Schema: CS421G68 SELECT C.NAME, COUNT(UNIQUE(CRID)) AS NUMBEROFRDATA...

22:28:51 SUCCESS SELECT 0.029 0.018 4 Success: 4 rows set fetched

22:28:51 FINISHED 0.029 0.018 4 Success: 4 rows set fetched

### 5.3 Who is worst at experiments?

Rank active researchers (name) by the number of raw data they produced that is associated with a dashboard containing focus > 3 (this needs to be normalized by the total number of experiments they have done):

```

20 SELECT M.NAME, M.NUMBEROFRDATA, N.NUMBEROFRDATA, ROUND(CAST(M.NUMBEROFRDATA AS DOUBLE)/N.NUMBEROFRDATA,3)
21 FROM
22 (SELECT Y.NAME, COUNT(X.CRID) AS NUMBEROFRDATA
23 FROM
24 (SELECT UNIQUE(D.CRID) AS CRID FROM CLICDASHBOARD D, CLICRAWDATA R
25 WHERE D.CRID = R.CRID AND D.FOCUS > 3) X, LABMEMBER Y, EXPERIMENT Z
26 WHERE Y.STATUS = 'active' AND Z.CRID = X.CRID AND Y.NAME = Z.NAME
27 GROUP BY Y.NAME
28 ORDER BY NUMBEROFRDATA),
29 (SELECT Y.NAME, COUNT(X.CRID) AS NUMBEROFRDATA
30 FROM
31 (SELECT UNIQUE(D.CRID) AS CRID FROM CLICDASHBOARD D, CLICRAWDATA R
32 WHERE D.CRID = R.CRID ) X, LABMEMBER Y, EXPERIMENT Z
33 WHERE Y.STATUS = 'active' AND Z.CRID = X.CRID AND Y.NAME = Z.NAME
34 GROUP BY Y.NAME
35 ORDER BY NUMBEROFRDATA)
36 WHERE M.NAME=N.NAME
37

```

22:45 [520] [INS]

Log 1: CLICDASHBOARD [6] X

SQL/Command

NAME	NUMBEROFRDATA	NUMBEROFRDATA	ROUND(CAST(NUMBEROFRDATA AS DOUBLE)/NUMBEROFRDATA,3)
Kim Metera	18	28	0.643
Sabrina Leslie	22	36	0.611
Francis Stabile	24	36	0.667
Radin Tahvildari	24	34	0.706
Daniel Berard	28	35	0.8
Cindy Shaheen	30	38	0.789

Time Status Command Exec Fetch Rows Message SQL/Command

13:23:35 STARTED Executing for: 'cs421g' [DB2 LUW], Schema: CS421G68 SELECT M.NAME, M.NUMBEROFRDATA, N.NUMBEROFRDATA, ROUND(CAST(M.NUMBEROFRDATA AS DOUBLE)/N.NUMBEROFRDATA,3)...

13:23:35 SUCCESS SELECT 0.019 0.002 6 Success: 6 rows set fetched

13:23:35 FINISHED 0.019 0.002 6 Success: 6 rows set fetched

### 5.4 Get me all my good data!! Meeting with PI tomorrow!!

List the paths of all DNA unwinding data whose binding analysis results have analysis type "Decent" or "Accurate":

```

23
24
25 SELECT PATHS
26 FROM DNAUNWINDINGDATA D, BINDINGANALYSISRESULTS B
27 WHERE D.CRID=B.CRID AND (B.ANALYSISTYPE IN ('decent','accurate'));
28
29
30

```

34:1 [798] [INS]

Log 1: DNAUNWINDINGDATA [140] X

SQL/Command

PATHS
C:\BindingAnalysisResults\4
C:\BindingAnalysisResults\5
C:\BindingAnalysisResults\7
C:\BindingAnalysisResults\8
C:\BindingAnalysisResults\10
C:\BindingAnalysisResults\11
C:\BindingAnalysisResults\12
C:\BindingAnalysisResults\13
C:\BindingAnalysisResults\14
C:\BindingAnalysisResults\15
C:\BindingAnalysisResults\16
C:\BindingAnalysisResults\17
C:\BindingAnalysisResults\18
C:\BindingAnalysisResults\20
C:\BindingAnalysisResults\21
C:\BindingAnalysisResults\22
C:\BindingAnalysisResults\24
C:\BindingAnalysisResults\25
C:\BindingAnalysisResults\26
C:\BindingAnalysisResults\29
C:\BindingAnalysisResults\31
C:\BindingAnalysisResults\32
C:\BindingAnalysisResults\33
C:\BindingAnalysisResults\36

Time Status Command Exec Fetch Rows Message SQL/Command

22:30:51 STARTED Executing for: 'cs421g' [DB2 LUW], Schema: CS421G68 SELECT PATHS ...

22:30:51 SUCCESS SELECT 0.024 0.020 154 Success: 154 rows set fetched

22:30:51 FINISHED 0.024 0.020 154 Success: 154 rows set fetched

## 5.5 Getting ready to publish!!

List the paths of all protocols that satisfy:

- they are used in generating nano-particle data with nano-particle analysis results, which produce average radius higher than 40;
- they are used in the experiments for more than 2 times (try number in experiment is bigger than or equal to 2)

```

29
30 SELECT DISTINCT P.PATHS
31 FROM PROTOCOL_P, EXPERIMENT E, NANOPARTICLEANALYSISRESULTS N
32 WHERE N.CRID = E.CRID AND N.AVGADIUS > 40 AND E.PROTOCOLID=P.PROTOCOLID AND E.TRYNUMBER >= 2;
33
34
35
27:67 [791] [INS]
Log 1: PROTOCOL [50] X

```

PATHS
1 C:\Protocol\102.pdf
2 C:\Protocol\104.pdf
3 C:\Protocol\105.pdf
4 C:\Protocol\111.pdf
5 C:\Protocol\119.pdf
6 C:\Protocol\125.pdf
7 C:\Protocol\128.pdf
8 C:\Protocol\130.pdf
9 C:\Protocol\137.pdf
10 C:\Protocol\142.pdf
11 C:\Protocol\143.pdf
12 C:\Protocol\145.pdf
13 C:\Protocol\15.pdf
14 C:\Protocol\152.pdf
15 C:\Protocol\154.pdf
16 C:\Protocol\169.pdf
17 C:\Protocol\17.pdf
18 C:\Protocol\18.pdf
19 C:\Protocol\181.pdf
20 C:\Protocol\188.pdf
21 C:\Protocol\190.pdf
22 C:\Protocol\193.pdf
23 C:\Protocol\196.pdf
24 C:\Protocol\197.pdf
25 C:\Protocol\2.pdf

Format: <Select a Cell>

Time	Status	Command	Exec	Fetch	Rows	Message	SQL/Command
22:31:44	STARTED					Executing for: 'c\$421g' [DB2 LUW], Schema: CS421G68	
22:31:45	SUCCESS	SELECT	0.048	0.018	50	Result set fetched	SELECT DISTINCT P.PATHS...
22:31:45	FINISHED		0.048	0.018	50	Success: 1	

## 5.6 The BIG Picture!

Researchers supervising undergraduates need to keep track of the progress of their assistants. For every active assistant, supervisors would like to see the average metrics regarding data dashboard in the last year. In addition, they also want to know how many datasets assistants have finished analyzing in order to compare the performance.

```

120 SELECT RNAME, AVGFOCUS, AVGSNR, AVGCONT, NBCRID
121 FROM
122 (SELECT *
123 FROM SUPERVISE S
124 JOIN CLICDASHBOARD D
125 (SELECT NAME, AVGFOCUS) AS AVGFOCUS, AVG(SNR) AS AVGSNR, AVG(CONTAMINANTS) AS AVGCONT, COUNT(DISTINCT CRID) AS NBCRID
126 FROM CLICDASHBOARD
127 WHERE YEAR(DATES)-YEAR(CURRENT_DATE) < 2
128 GROUP BY NAME)
129 ON NAME=S.RNAME, LABMEMBER L
130 WHERE A.RNAME=L.NAME AND L.STATUS='active'
131 ORDER BY RNAME, RNAME;
132
119:1 [2917] [INS]
Log 1: SUPERVISE [4] X

```

RNAME	RNAME	AVGFOCUS	AVGSNR	AVGCONT	NBCRID
1 Cindy Shaheen	Stone Chen	3.1614285714285715	5.633214285714286	2272	46
2 Cindy Shaheen	Zhi Zhang	2.5334090909090907	5.752272727272728	2468	35
3 Radin Tahvildari	Wendy Ji	3.0553061224489797	5.5916326530612235	2526	46
4 Radin Tahvildari	Zach Friedenberger	3.108653846153845	5.386923076923078	2506	46

Time Status Command Exec Fetch Rows Message SQL/Command

23:03:54 STARTED Executing for: 'c\$421g' [DB2 LUW], Schema: CS421G68

23:03:54 SUCCESS SELECT 0.028 0.024 4 Result set fetched

23:03:54 FINISHED 0.028 0.024 4 Success: 1

SELECT RNAME, RNAME, AVGFOCUS, AVGSNR, AVGCONT, NBCRID...

## 6 Data Modification

### 6.1 Research assistant has graduated...

Find the research assistants named 'Stone Chen' and 'Wendy Ji' and modify their status to 'inactive'. For those research assistants who just left, update their end dates of supervision to "2019-05-01":

- Table LABMEMBER:

Previous Data:

```
112 SELECT * FROM LABMEMBER;
113
114
115
108:29 [2957] INS [ ]
```

*	NAME	CONTACTINFORMATION	STATUS
1	Wendy Ji	(514)-561-5587	active
2	Stone Chen	(514)-531-3322	active
3	Zhi Zhang	(514)-557-3232	active
4	Zach Friedenberger	(514)-324-1237	active
5	Gracia Gu	(514)-545-5123	active
6	Zhi Yue Zhang	(514)-512-5567	active
7	Cindy Shaheen	(514)-532-8521	active
8	Francis Stabile	(514)-274-4887	active
9	Daniel Berard	(514)-244-9817	active
10	Kim Metera	(514)-274-4921	active
11	Radin Tahvildari	(514)-284-1237	active
12	Romain Berti	(514)-224-9387	active
13	Sabrina Leslie	(514)-224-2837	active
14	Shane Scott	(514)-274-9287	inactive
15	Haoran Liao	(514)-271-1927	inactive
16	Yash Patel	(514)-292-1927	inactive

Data after modification:

```
38 UPDATE LABMEMBER
39 SET STATUS = CASE
40 WHEN NAME IN ('Stone Chen', 'Wendy Ji') THEN 'inactive'
41 ELSE STATUS
42 END;
43
44 SELECT * FROM LABMEMBER;
45
35:1 [977] INS [ ]
```

*	NAME	CONTACTINFORMATION	STATUS
1	Wendy Ji	(514)-561-5587	inactive
2	Stone Chen	(514)-531-3322	inactive
3	Zhi Zhang	(514)-557-3232	active
4	Zach Friedenberger	(514)-324-1237	active
5	Gracia Gu	(514)-545-5123	active
6	Zhi Yue Zhang	(514)-512-5567	active
7	Cindy Shaheen	(514)-532-8521	active
8	Francis Stabile	(514)-274-4887	active
9	Daniel Berard	(514)-244-9817	active
10	Kim Metera	(514)-274-4921	active
11	Radin Tahvildari	(514)-284-1237	active
12	Romain Berti	(514)-224-9387	active
13	Sabrina Leslie	(514)-224-2837	active

- Table SUPERVISE:

Previous Data:

```
61 SELECT * FROM SUPERVISE;
62
63
54:19 [1373] INS [ ]
```

*	STARTDATE	ENDDATE	RNAME	RANAME
1	2018-05-01	(null)	Radin Tahvildari	Wendy Ji
2	2018-05-01	(null)	Radin Tahvildari	Zach Friedenberger
3	2017-02-01	2018-06-01	Radin Tahvildari	Haoran Liao
4	2017-02-01	2018-06-01	Radin Tahvildari	Yash Patel
5	2018-05-01	(null)	Cindy Shaheen	Stone Chen
6	2018-05-01	(null)	Cindy Shaheen	Zhi Zhang
7	2018-05-01	(null)	Cindy Shaheen	Gracia Gu

Data after modification:

```
46 UPDATE BINDINGANALYSISRESULTS
47 SET ANALYSISTYPE = 'decent'
48 WHERE NAME='Zhi Zhang' AND ANALYSISTYPE = 'rough';
49
50 UPDATE SUPERVISE S
51 SET ENDDATE = '2019-05-01'
52 WHERE EXISTS
53 (SELECT L.NAME
54 FROM LABMEMBER L
55 WHERE L.STATUS = 'inactive' AND L.NAME=S.RNAME AND S.ENDDATE IS NULL);
56
57 SELECT * FROM SUPERVISE;
58
```

*	STARTDATE	ENDDATE	RNAME	RANAME
1	2018-05-01	2019-05-01	Radin Tahvildari	Wendy Ji
2	2018-05-01	(null)	Radin Tahvildari	Zach Friedenberger
3	2017-02-01	2018-06-01	Radin Tahvildari	Haoran Liao
4	2017-02-01	2018-06-01	Radin Tahvildari	Yash Patel
5	2018-05-01	2019-05-01	Cindy Shaheen	Stone Chen
6	2018-05-01	(null)	Cindy Shaheen	Zhi Zhang
7	2018-05-01	(null)	Cindy Shaheen	Gracia Gu

## 6.2 Newbie to the lab!!

Fortunately, someone is replacing those who left. Insert into research assistant and lab member "Raffle Zhu" with phone number (514)-888-8888 with 'active' status and supervision start date "2019-04-01" by "Zhi Yue Zhang":

Database Connection      Sticky Database      Schema

```

1 --SELECT End'Raffle' AS NEW ,NULL AS N FROM SUPERVISE WHERE RENAME='Stone Chen';
2
3 INSERT INTO LABMEMBER
4 VALUES ('Raffle Zhu','(514)-888-8888','active');
5
6 INSERT INTO RESEARCHASSISTANT
7 VALUES ('Raffle Zhu');
8
9 INSERT INTO SUPERVISE (STARTDATE, RNAME, RENAME)
10 VALUES ('2019-04-01', (SELECT RNAME FROM SUPERVISE
11 WHERE RENAME='Stone Chen'), 'Raffle Zhu');
12
13

```

1:4 [q] [INS]

Log

Time	Status	Command	Exec	Fetch	Rows	Message
12:29:58	STARTED					Executing for: 'cs421' [DB2 LUW], Schema: CS421G68
→ 12:29:58	SUCCESS	INSERT	0.040		1	OK
▲ 12:29:59	FINISHED		0.040	0	1	Success: 1

Previous Data:

Table: LABMEMBER

cs421/CS421G68/TABLE/LABMEMBER

Info Columns Data Row Count Primary Key

In DbVisualizer Pro this feature includes data editing functionality and management.

*	NAME	CONTACTINFORMATION	STATUS
1	Wendy Ji	(514)-561-5587	inactive
2	Stone Chen	(514)-531-3322	inactive
3	Zhi Zhang	(514)-557-3232	active
4	Zach Friedenberger	(514)-324-1237	active
5	Gracia Gu	(514)-545-5123	active
6	Zhi Yue Zhang	(514)-512-5567	active
7	Cindy Shaheen	(514)-532-8521	active
8	Francis Stabile	(514)-274-4887	active
9	Daniel Berard	(514)-244-9817	active
10	Kim Metera	(514)-274-4921	active
11	Radin Tahvildari	(514)-284-1237	active
12	Romain Berti	(514)-224-9387	active
13	Sabrina Leslie	(514)-224-2837	active
14	Shane Scott	(514)-274-9287	inactive
15	Haoran Liao	(514)-271-1927	inactive
16	Yash Patel	(514)-292-1927	inactive

Data after modification:

Table: LABMEMBER

cs421/CS421G68/TABLE/LABMEMBER

Info Columns Data Row Count Primary Key

In DbVisualizer Pro this feature includes data editing functionality and management.

*	NAME	CONTACTINFORMATION	STATUS
1	Wendy Ji	(514)-561-5587	inactive
2	Stone Chen	(514)-531-3322	inactive
3	Zhi Zhang	(514)-557-3232	active
4	Zach Friedenberger	(514)-324-1237	active
5	Gracia Gu	(514)-545-5123	active
6	Zhi Yue Zhang	(514)-512-5567	active
7	Cindy Shaheen	(514)-532-8521	active
8	Francis Stabile	(514)-274-4887	active
9	Daniel Berard	(514)-244-9817	active
10	Kim Metera	(514)-274-4921	active
11	Radin Tahvildari	(514)-284-1237	active
12	Romain Berti	(514)-224-9387	active
13	Sabrina Leslie	(514)-224-2837	active
14	Shane Scott	(514)-274-9287	inactive
15	Haoran Liao	(514)-271-1927	inactive
16	Yash Patel	(514)-292-1927	inactive
17	Raffle Zhu	(514)-888-8888	active

Previous Data:

Table: SUPERVISE

cs421/CS421G68/TABLE/SUPERVISE

Info Columns Data Row Count Primary Key Index

In DbVisualizer Pro this feature includes data editing functionality and management.

*	STARTDATE	ENDDATE	RNAME	RANAME
1	2018-05-01	(null)	Radin Tahvildari	Wendy Ji
2	2018-05-01	(null)	Radin Tahvildari	Zach Friedenberger
3	2017-02-01	2018-06-01	Radin Tahvildari	Haoran Liao
4	2017-02-01	2018-06-01	Radin Tahvildari	Yash Patel
5	2018-05-01	(null)	Cindy Shaheen	Stone Chen
6	2018-05-01	(null)	Cindy Shaheen	Zhi Zhang
7	2018-05-01	(null)	Cindy Shaheen	Gracia Gu

Data after modification:

Table: SUPERVISE

cs421/CS421G68/TABLE/SUPERVISE

Info Columns Data Row Count Primary Key Index

In DbVisualizer Pro this feature includes data editing functionality and management.

*	STARTDATE	ENDDATE	RNAME	RANAME
1	2018-05-01	(null)	Radin Tahvildari	Wendy Ji
2	2018-05-01	(null)	Radin Tahvildari	Zach Friedenberger
3	2017-02-01	2018-06-01	Radin Tahvildari	Haoran Liao
4	2017-02-01	2018-06-01	Radin Tahvildari	Yash Patel
5	2018-05-01	(null)	Cindy Shaheen	Stone Chen
6	2018-05-01	(null)	Cindy Shaheen	Zhi Zhang
7	2018-05-01	(null)	Cindy Shaheen	Gracia Gu
8	2019-04-01	(null)	Cindy Shaheen	Raffle Zhu

Previous Data:

Table: RESEARCHASSISTANT					
cs421/CS421G68/TABLE/RESEARCHASSISTANT					
		Info	Columns	Data	Row Count
<b>i</b> In DbVisualizer Pro this feature includes data					
<b>*</b>	NAME				
1	Gracia Gu				
2	Haoran Liao				
3	Stone Chen				
4	Wendy Ji				
5	Yash Patel				
6	Zach Friedenberger				
7	Zhi Zhang				

Data after modification:

Table: RESEARCHASSISTANT					
cs421/CS421G68/TABLE/RESEARCHASSISTANT					
		Info	Columns	Data	Row Count
<b>i</b> In DbVisualizer Pro this feature includes data					
<b>*</b>	NAME				
1	Gracia Gu				
2	Haoran Liao				
3	Raffle Zhu				
4	Stone Chen				
5	Wendy Ji				
6	Yash Patel				
7	Zach Friedenberger				
8	Zhi Zhang				

### 6.3 Correct a mistake!!

Find all DNA binding analysis results done by "Zhi Zhang" in 2019 with analysis type "rough". Correct the types of those analysis to "decent":

```

72 DELETE FROM CLICASHBOARD
73 WHERE FOCUS < 2 OR (YEAR(CURRENT_DATE)-YEAR(DATES)) > 1;
74
75
76 SELECT * FROM CLICASHBOARD;
66:31 [1579] INS
Log
Time Status Command Exec Fetch Rows Message SQL/Command
12:56:59 STARTED Executing for: `cs421g` [DB2 LUW], Schema: CS421G68
12:56:59 SUCCESS DELETE 0.012 57 OK DELETE FROM CLICASHBOARD...
12:56:59 FINISHED 0.012 0 57 Success: 1

```

Previous Data:

76 SELECT * FROM CLICASHBOARD;								
78:20 [1915] INS								
Log 1: CLICASHBOARD [201] x								
*	DID	DATES	PATHS	FOCUS	SNR	CONTAMINANTS	CRID	NAME
189	189	2018-08-28	C:\CLICDashboard\189	4.64	4.09	2900	188	Zhi Zhang
190	190	2018-06-07	C:\CLICDashboard\190	4.39	2.99	459	67	Stone Chen
191	191	2019-01-29	C:\CLICDashboard\191	2.36	4.97	2307	27	Wendy Ji
192	192	2019-02-20	C:\CLICDashboard\192	2.84	7.62	2127	189	Stone Chen
193	193	2018-12-20	C:\CLICDashboard\193	3.09	5.85	4671	47	Stone Chen
194	194	2018-07-05	C:\CLICDashboard\194	1.73	6.64	4948	105	Zhi Zhang
195	195	2018-10-22	C:\CLICDashboard\195	3.84	8.71	2359	194	Wendy Ji
196	196	2018-07-14	C:\CLICDashboard\196	4.8	5.92	4848	50	Stone Chen
197	197	2018-07-06	C:\CLICDashboard\197	4.75	8.12	4942	87	Stone Chen
198	198	2019-01-16	C:\CLICDashboard\198	1.07	7.27	2758	170	Zhi Zhang
199	199	2019-01-12	C:\CLICDashboard\199	3.87	8.73	3283	114	Zhi Zhang
200	200	2018-07-21	C:\CLICDashboard\200	2.56	4.98	3778	159	Zach Friedenberger
201	202	2016-03-03	C:\CLICDashboard\202	4.0	3.14	1000	115	Stone Chen

Data after modify:

72 DELETE FROM CLICASHBOARD								
73 WHERE FOCUS < 2 OR (YEAR(CURRENT_DATE)-YEAR(DATES)) > 1;								
74								
75								
*	DID	DATES	PATHS	FOCUS	SNR	CONTAMINANTS	CRID	NAME
132	185	2018-05-26	C:\CLICDashboard\185	2.39	8.36	2333	125	Wendy Ji
133	187	2019-01-31	C:\CLICDashboard\187	2.95	3.4	3044	157	Zach Friedenberger
134	188	2019-02-09	C:\CLICDashboard\188	3.41	8.12	4003	153	Wendy Ji
135	189	2018-08-28	C:\CLICDashboard\189	4.64	4.09	2900	188	Zhi Zhang
136	190	2018-06-07	C:\CLICDashboard\190	4.39	2.99	459	67	Stone Chen
137	191	2019-01-29	C:\CLICDashboard\191	2.36	4.97	2307	27	Wendy Ji
138	192	2019-02-26	C:\CLICDashboard\192	2.84	7.62	2127	189	Stone Chen
139	193	2018-12-20	C:\CLICDashboard\193	3.09	5.85	4671	47	Stone Chen
140	195	2018-10-22	C:\CLICDashboard\195	3.84	8.71	2359	194	Wendy Ji
141	196	2018-07-14	C:\CLICDashboard\196	4.8	5.92	4848	50	Stone Chen
142	197	2018-07-06	C:\CLICDashboard\197	4.75	8.12	4942	87	Stone Chen
143	199	2019-01-12	C:\CLICDashboard\199	3.87	8.73	3283	114	Zhi Zhang
144	200	2018-07-21	C:\CLICDashboard\200	2.56	4.98	3778	159	Zach Friedenberger

### 6.4 Clean up some trash data!!

We delete those useless data (satisfy both of the following conditions):

- Data with serious errors (the focus is less than 2);

- Data generated long ago (has been generated for more than 1 years)

Previous Data:

```

70
71
72
73
74 SELECT * FROM CLICDASHBOARD;
75
76
77
78
79
80
81
82
83
74:29 [1923] [INS]
Log 1: CLICDASHBOARD [202] x
Format: <Select a Cell>

```

*	DID	DATES	PATHS	FOCUS	SNR	CONTAMINANTS	CRID	NAME
183	183	2018-01-08	C:\CLICDashboard\183	3.09	9.14	1912	81	Wendy Ji
184	184	2019-02-10	C:\CLICDashboard\184	1.43	8.61	1225	63	Wendy Ji
185	185	2018-05-26	C:\CLICDashboard\185	2.39	8.36	2333	125	Wendy Ji
186	186	2018-06-03	C:\CLICDashboard\186	1.71	2.22	1105	63	Zhi Zhang
187	187	2019-01-31	C:\CLICDashboard\187	2.95	3.4	3044	197	Zach Friedenberger
188	188	2019-02-09	C:\CLICDashboard\188	3.41	8.12	4003	153	Wendy Ji
189	189	2018-08-28	C:\CLICDashboard\189	4.64	4.09	2900	188	Zhi Zhang
190	190	2018-06-07	C:\CLICDashboard\190	4.39	2.99	459	67	Stone Chen
191	191	2019-01-29	C:\CLICDashboard\191	2.36	4.97	2307	27	Wendy Ji
192	192	2019-02-26	C:\CLICDashboard\192	2.84	7.62	2127	189	Stone Chen
193	193	2018-12-20	C:\CLICDashboard\193	3.09	5.85	4671	47	Stone Chen
194	194	2018-07-05	C:\CLICDashboard\194	1.73	6.64	4948	105	Zhi Zhang
195	195	2018-10-22	C:\CLICDashboard\195	3.84	8.71	2359	194	Wendy Ji
196	196	2018-07-14	C:\CLICDashboard\196	4.8	5.92	4848	50	Stone Chen
197	197	2018-07-06	C:\CLICDashboard\197	4.75	8.12	4942	87	Stone Chen
198	198	2019-01-16	C:\CLICDashboard\198	1.07	7.27	2758	170	Zhi Zhang
199	199	2019-01-12	C:\CLICDashboard\199	3.87	8.73	3283	114	Zhi Zhang
200	200	2018-07-21	C:\CLICDashboard\200	2.56	4.98	3778	159	Zach Friedenberger
201	201	2016-03-03	C:\CLICDashboard\201	1.0	0.67	5000	115	Stone Chen
202	202	2016-03-03	C:\CLICDashboard\202	4.0	3.14	1000	115	Stone Chen

Data after modification:

```

64
65
66
67
68 DELETE FROM CLICDASHBOARD
69 WHERE CONTAMINANTS/FOCUS > 4500 AND (YEAR(CURRENT_DATE)-YEAR(DATES)) > 2;
70
71
72 SELECT * FROM CLICDASHBOARD;
73
74
75
76
73:1 [1922] [INS]
Log 1: CLICDASHBOARD [201] x
Format: <Select a Cell>

```

*	DID	DATES	PATHS	FOCUS	SNR	CONTAMINANTS	CRID	NAME
182	182	2018-12-16	C:\CLICDashboard\182	1.84	3.76	822	141	Zhi Zhang
183	183	2019-01-08	C:\CLICDashboard\183	3.09	9.14	1912	81	Wendy Ji
184	184	2019-02-10	C:\CLICDashboard\184	1.43	8.61	1225	63	Wendy Ji
185	185	2018-05-26	C:\CLICDashboard\185	2.39	8.36	2333	125	Wendy Ji
186	186	2018-06-03	C:\CLICDashboard\186	1.71	2.22	1105	63	Zhi Zhang
187	187	2019-01-31	C:\CLICDashboard\187	2.95	3.4	3044	197	Zach Friedenberger
188	188	2019-02-09	C:\CLICDashboard\188	3.41	8.12	4003	153	Wendy Ji
189	189	2018-08-28	C:\CLICDashboard\189	4.64	4.09	2900	188	Zhi Zhang
190	190	2018-06-07	C:\CLICDashboard\190	4.39	2.99	459	67	Stone Chen
191	191	2019-01-29	C:\CLICDashboard\191	2.36	4.97	2307	27	Wendy Ji
192	192	2019-02-26	C:\CLICDashboard\192	2.84	7.62	2127	189	Stone Chen
193	193	2018-12-20	C:\CLICDashboard\193	3.09	5.85	4671	47	Stone Chen
194	194	2018-07-05	C:\CLICDashboard\194	1.73	6.64	4948	105	Zhi Zhang
195	195	2018-10-22	C:\CLICDashboard\195	3.84	8.71	2359	194	Wendy Ji
196	196	2018-07-14	C:\CLICDashboard\196	4.8	5.92	4848	50	Stone Chen
197	197	2018-07-06	C:\CLICDashboard\197	4.75	8.12	4942	87	Stone Chen
198	198	2019-01-16	C:\CLICDashboard\198	1.07	7.27	2758	170	Zhi Zhang
199	199	2019-01-12	C:\CLICDashboard\199	3.87	8.73	3283	114	Zhi Zhang
200	200	2018-07-21	C:\CLICDashboard\200	2.56	4.98	3778	159	Zach Friedenberger
201	201	2016-03-03	C:\CLICDashboard\201	1.0	0.67	5000	115	Stone Chen
202	202	2016-03-03	C:\CLICDashboard\202	4.0	3.14	1000	115	Stone Chen

## 7 Views

### 7.1 You are FIRED!!

Of course this does not happen and this is just a simulation. Rank all active research assistants who started in 2018 (check this with supervision start date) by the number of data they have analyzed (including both types of data):

```

73
74 DROP VIEW SOMENAME;
75
76 CREATE VIEW SOMENAME (NAME) AS
77 SELECT B.NAME
78 FROM BINDINGANALYSISRESULTS B, SUPERVISE S, LABMEMBER L
79 WHERE B.NAME = S.RNAME AND YEAR(S.STARTDATE) >= 2018 AND B.NAME = L.NAME AND L.STATUS = 'active'
80 UNION ALL
81 SELECT N.NAME
82 FROM NANOPARTICLEANALYSISRESULTS N, SUPERVISE S, LABMEMBER L
83 WHERE N.NAME = S.RNAME AND YEAR(S.STARTDATE) >= 2018 AND N.NAME = L.NAME AND L.STATUS = 'active';
84
85 SELECT * FROM SOMENAME;
86 SELECT NAME,COUNT(*) AS NUMBERDA
87 FROM SOMENAME
88 GROUP BY NAME
89 ORDER BY NUMBERDA;
90
91
92
93
94
95
96
97
86:33 [2385] [INS]
Log 1: SOMENAME [4] x

```

*	NAME	NUMBERDA
1	Wendy Ji	93
2	Zach Friedenberger	93
3	Zhi Zhang	97
4	Stone Chen	117

Time	Status	Command	Exec	Fetch	Rows	Message	SQI/Command
▼ 22:33:19	■ STARTED	Executing for: 'c5421g' [DB2 LUW], Schema: C5421G68					
→ 22:33:20	✓ SUCCESS	DROP	0.026		0	OK. No rows were affected	DROP VIEW SOMENAME
→ 22:33:20	✓ SUCCESS	CREATE	0.026		0	OK. No rows were affected	CREATE VIEW SOMENAME (NAME) AS...
→ 22:33:20	✓ SUCCESS	SELECT	0.044	0.022	400	Result set fetched	SELECT * FROM SOMENAME
→ 22:33:20	✓ SUCCESS	SELECT	0.124	0.020	4	Result set fetched	SELECT NAME,COUNT(*) AS NUMBERDA...
▲ 22:33:20	■ FINISHED		0.220	0.042	404	Success: 4	

The View SOMENAME is not updatable since the attribute list of the view does not include the primary key, which makes the update unable to locate the entry accurately.

## 7.2 Time to check the equipment...

Laboratory equipment needs to be checked and maintained frequently to ensure the smooth running of the experiment and the accuracy of its result, so we will count the number of times a microscope has been used. One will check and maintain this equipment when the number of uses exceeds 1000 and reset the number of uses to 0 after maintenance.

```

91 DROP VIEW MICROSCOPECONDITIONS;
92
93 CREATE VIEW MICROSCOPECONDITIONS (MICROID,NUMBERUSED) AS
94 SELECT MICROSOPED,COUNT(*) AS NUMBERUSED
95 FROM
96 (SELECT E.MICROSOPED
97 FROM EXPERIMENT E, MICROSCOPE M
98 WHERE E.MICROSOPED = M.MICROSOPED AND M.STATUS = 'active')T
99 GROUP BY MICROSOPED
100 ;
101
102 SELECT * FROM MICROSCOPECONDITIONS;
103
104
105 UPDATE MICROSCOPE M
106 SET M.STATUS = 'inactive'
107 WHERE EXISTS
108 (SELECT *
109 FROM MICROSCOPECONDITIONS C
110 WHERE C.NUMBERUSED >= 1000);
111

```

*	MICROID	NUMBERUSED
1	Nikon01	195
2	Nikon02	187

The View MICROSCOPECONDITIONS is not updatable it contains the GROUP BY clause:

## 7.3

The View is updatable when it satisfies all the following conditions:

- The view is defined based on only one table and includes its PRIMARY KEY upon creation.
- The view should not have any field made out of aggregate and arithmetic functions.
- The view must not have any GROUP BY, HAVING, SUBQUERIES or DISTINCT clause in its definition.
- If the view you want to update is based upon another view, the later should be updatable.
- Any of the selected output fields (of the view) must not use constants, strings or value expressions.
- (DB2) The view is deletable and not view only.
- (DB2) UNION ALL must have matching types for each field.

## 8 CHECK Constraint

### 8.1 Date check

Add the check constraint to the table BINDINGANALYSISRESULTS, to make sure that all the data stored in the database are after Jan, 2018.

```
130 ALTER TABLE BINDINGANALYSISRESULTS
131 ADD CONSTRAINT DATE_check CHECK (YEAR(DATES) >= 2018);
132
133 ALTER TABLE BINDINGANALYSISRESULTS
134 DROP CONSTRAINT DATE_check;
135
136
```

126:35 [3410] [INS]

Time	Status	Command	Exec	Fetch	Rows	Message	SQL/Command
22:56:32	STARTED					Executing for: 'cs421g' [DB2 LUW], Schema: CS421G68	
22:56:32	SUCCESS	ALTER	0.035	0	0	No rows were affected	ALTER TABLE BINDINGANALYSISRESULTS ...
22:56:32	SUCCESS	ALTER	0.054	0	0	No rows were affected	ALTER TABLE BINDINGANALYSISRESULTS ...
22:56:32	FINISHED		0.089	0	0	Success: 2	

We tried to insert a tuple whose date is earlier than 2018, leading to failure of the sql execution.

```
155 INSERT INTO BINDINGANALYSISRESULTS
156 VALUES ('201, '2010-09-02', 'C:\BindingAnalysisResults\201', 'accurate', 83, 'Stone_Chen');
157
158 ALTER TABLE BINDINGANALYSISRESULTS
159 DROP CONSTRAINT DATE_CHECK;
160
161
162
```

143:20 [3547] [INS]

Time	Status	Command	Exec	Fetch	Rows	Message	SQL/Command
09:55:42	STARTED					Executing for: 'cs421g' [DB2 LUW], Schema: CS421G68	
09:55:42	FAILED	INSERT	0.112	0	0	0 Code: -545, SQL State: 23513] The requested operation is not allowed because a row does not satisfy the check constraint 'CS421G68.BINDINGANALYSISRESULTS.DATE_CHECK'.	INSERT INTO BINDINGANALYSISRESULTS VALUES ('201, '2010-09-02', 'C:\BindingAnalysisResults\201', 'accurate', 83, 'Stone_Chen')
09:55:42	FINISHED						

### 8.2 Number check

Add the check constraint to the contact information (the phone number) in the table LABMEMBER to make sure the phone number is valid.

```
136 ALTER TABLE LABMEMBER
137 ADD CONSTRAINT NUMBER_check CHECK (CONTACTINFORMATION LIKE '(\%)-\d{3}-\d{4}');
138
139 ALTER TABLE LABMEMBER
140 DROP CONSTRAINT NUMBER_check;
141
142
```

130:26 [3471] [INS]

Time	Status	Command	Exec	Fetch	Rows	Message	SQL/Command
23:01:23	STARTED					Executing for: 'cs421g' [DB2 LUW], Schema: CS421G68	
23:01:23	SUCCESS	ALTER	0.036	0	0	OK. No rows were affected	ALTER TABLE LABMEMBER ...
23:01:23	SUCCESS	ALTER	0.030	0	0	OK. No rows were affected	ALTER TABLE LABMEMBER ...
23:01:23	FINISHED		0.066	0	0	Success: 2	

We tried to insert a tuple which violated our phone number constraint, leading to an failure of the sql command.

```
165 ALTER TABLE LABMEMBER
166 ADD CONSTRAINT NUMBER_check CHECK (CONTACTINFORMATION LIKE '(\%)-\d{3}-\d{4}');
167
168 --TEST CASE
169 INSERT INTO LABMEMBER
170 VALUES ('Mark Zhu', '(514)8888888', 'active');
171
172 --TEST
173 ALTER TABLE LABMEMBER
174 DROP CONSTRAINT NUMBER_check;
175
```

152:23 [3550] [INS]

Time	Status	Command	Exec	Fetch	Rows	Message	SQL/Command
09:52:27	STARTED					Executing for: 'cs421g' [DB2 LUW], Schema: CS421G68	
09:52:27	FAILED	INSERT	0.310	0	0	0 Code: -545, SQL State: 23513] The requested operation is not allowed because a row does not satisfy the check constraint 'CS421G68.LABMEMBER.NUMBER_CHECK'.	INSERT INTO LABMEMBER VALUES ('Mark Zhu', '(514)8888888', 'active')
09:52:27	FINISHED						

## 9 Creative Points

We have implemented or performed the followings to our database:

- Automatic Data Generation:** We wrote a python script to automatically construct SQL insert statements with randomly generated data. Please refer to autoDataGen.py and autoDataGen.sql.
- Some Real Data:** We incorporated some real data in our database, including lab members and their corresponding status, as well as microscope information. In terms of other automatically generated data, they also follow the format of real data.
- Complex Analytical Query:** All queries in section 5 are real life questions one may ask in the lab and some of them involve certain complexities.

# CLiC Microscopic Data Organization Database P3(Group 68)

Stone Chen, Yunwen Ji, Wenquan Li, Yuki Zhang

March 21, 2019

## 1 Stored Procedure

### 1.1 Dubious Results

We created a stored procedure to if the average contaminants generated from a data source is above an acceptable threshold. If above, the procedure would change the analysis result type conducted on this data from 'decent' to 'rough' since the data quality was questionable.

```
cs421=> SELECT * from bindinganalysisresults where crid=1;
      barid |   dates   |          paths          | analysistype | crid |      name
-----+-----+-----+-----+-----+-----+-----+
      183 | 2018-12-14 | C:\BindingAnalysisResults\183 | decent      |     1 | Zach Friedenberger
      102 | 2018-10-05 | C:\BindingAnalysisResults\102 | decent      |     1 | Wendy Ji
(2 rows)

cs421=> SELECT * FROM CLICDASHBOARD WHERE CRID=1;
      did |   dates   |          paths          | focus | snr | contaminants | crid |      name
-----+-----+-----+-----+-----+-----+-----+
       9 | 2018-07-24 | C:\CLiCDashBoard\9 |  4.01 | 1.4 |        773 |     1 | Wendy Ji
      31 | 2019-01-14 | C:\CLiCDashBoard\31 |  4.74 | 2.16 |      3849 |     1 | Zach Friedenberger
     147 | 2018-08-18 | C:\CLiCDashBoard\147 |  4.39 | 5.03 |      2231 |     1 | Zach Friedenberger
(3 rows)

cs421=> create or replace function downgrade (target integer, threshold numeric) returns void as $$%
cs421$> declare
cs421$>   c refcursor;
cs421$>   c_rec record;
cs421$>   num integer;
cs421$>   average numeric;
cs421$> begin
cs421$>   num := 0;
cs421$>   average := 0;
cs421$>   open c for select * from clicdashboard where crid=target;
cs421$>   loop
cs421$>     fetch c into c_rec;
cs421$>     exit when not found;
cs421$>     if c_rec.contaminants is not null then
cs421$>       num = num + 1;
cs421$>       average = (average*(num-1)+c_rec.contaminants)/num;
cs421$>     end if;
cs421$>   end loop;
cs421$>   close c;
cs421$>   if average > threshold and threshold > 0 then
cs421$>     update bindinganalysisresults set analysistype='rough' where crid=target and analysistype='decent';
cs421$>   end if;
cs421$> end;
cs421$> $$ language plpgsql;
CREATE FUNCTION
cs421=> select downgrade(1,2000);
      downgrade
-----
(1 row)

cs421=> select * from bindinganalysisresults where crid = 1;
      barid |   dates   |          paths          | analysistype | crid |      name
-----+-----+-----+-----+-----+-----+-----+
      183 | 2018-12-14 | C:\BindingAnalysisResults\183 | rough      |     1 | Zach Friedenberger
      102 | 2018-10-05 | C:\BindingAnalysisResults\102 | rough      |     1 | Wendy Ji
(2 rows)
```

As we can see, the created function changed the analysis type of the two results related to the problematic data source with crid=1 because the average contaminant was above 2000.

## 2 GUI Based Data Retriever

We chose to do a demo on Mar 27, 2019 1:50-2:00 P.M. Our application was designed to allow database user to conveniently retrieve information on the CLiC data and its relevant metrics. In particular, the application is GUI based, consisted of panels with quick actions, custom queries as well as a CLiC data filter. All errors are handled at the back-end and error messages received from the database are displayed on the GUI. For most of the user input, especially in the CLiC data filter, input validity is assessed real time with feedback displayed.

## 3 Index

### 3.1 Chronological is Logical

We created index on DATES of CLICDASHBOARD.

Previous Data:

```
Query Editor Query History
1 SELECT C.NAME, COUNT(CRID) AS NUMBEROFTDATA
2 FROM CLICDASHBOARD C, LABMEMBER L
3 WHERE C.DATES >= '2019-01-01' AND C.DATES <= '2019-02-28'
4 AND C.NAME = L.NAME AND L.STATUS='active'
5 GROUP BY C.NAME
6 ORDER BY NUMBEROFTDATA;
7
8
9
10
11
12
13
14
15
16
17
```

Data Output Explain Messages Notifications

Successfully run. Total query runtime: 539 msec.  
3 rows affected.

After:

```
Query Editor Query History
1 SELECT C.NAME, COUNT(CRID) AS NUMBEROFTDATA
2 FROM CLICDASHBOARD C, LABMEMBER L
3 WHERE C.DATES >= '2019-01-01' AND C.DATES <= '2019-02-28'
4 AND C.NAME = L.NAME AND L.STATUS='active'
5 GROUP BY C.NAME
6 ORDER BY NUMBEROFTDATA;
7
8
9 CREATE INDEX date_idx ON CLICDASHBOARD (DATES);
10 DROP INDEX date_idx;
11
12
13
14
15
16
17
```

Data Output Explain Messages Notifications

Successfully run. Total query runtime: 88 msec.  
3 rows affected.

The results without the indexes were returned in 539 msec, while the results with the indexes were took 88 msec.

It is useful if we want to find the dates in a certain range and it should shorten the time required for the queries.

This is a considerable difference, because entries can help quickly locate data without having to search every row in a database table every time a database table is accessed, so our query is speeding up.

### 3.2 We Value Personal Achievements

We create index on NAME of NANOPARTICLEANALYSISRESULTS and BINDINGANALYSISRESULTS.

Previous Data:

```
Query Editor Query History
19
20
21 SELECT PATHS
22 FROM NANOPARTICLEANALYSISRESULTS
23 WHERE NAME = 'Wendy_31';
24 UNION
25 SELECT PATHS
26 FROM BINDINGANALYSISRESULTS
27 WHERE NAME = 'Wendy_31';
28
29
30
31
32
33
34
35
```

Data Output Explain Messages Notifications

Successfully run. Total query runtime: 393 msec.  
93 rows affected.

After:

```
Query Editor  Query History
18
19
20
21 SELECT PATHS
22 FROM NANOPARTICLEANALYSISRESULTS
23 WHERE NAME = 'Wendy Ji'
24 UNION
25 SELECT PATHS
26 FROM BINDINGANALYSISRESULTS
27 WHERE NAME = 'Wendy Ji';
28
29 CREATE INDEX info_idx ON NANOPARTICLEANALYSISRESULTS (NAME);
30 CREATE INDEX info_idx ON BINDINGANALYSISRESULTS (NAME);
31 DROP INDEX info_idx;
32 DROP INDEX info_idx;
33
34
Data Output Explain Messages Notifications
Successfully run. Total query runtime: 159 msec.
93 rows affected.
```


The results without the indexes were 393 msec, while the results with the indexes were took 159 msec. It is useful because we want to find the paths of a certain name which is similar to locate data, so indexes can also help us to reduce the time in having to search every row in tables. Therefore our query also speed up.

## 4 Visualization

### 4.1 Who is the hardest worker?

Rank the total number of results analyzed by the four research assistant who joined the lab at the same time, where results were counted from both nanoparticle analysis and binding analysis:


```
|cs421=> \COPY (SELECT NAME, COUNT(*) FROM (SELECT NAME FROM NANOPARTICLEANALYSISRESULTS UNION ALL SELECT NAME FROM BINDINGANALYSISRESULTS)T GROUP BY NAME)TO results.csv WITH CSV
COPY 4 _
```



### 4.2 Focus Measurement as Experiment Quality Metric

Get the focus measurements of Cindy Shaheen's experiment to find out whether the focus in her experiment is stable or not, indicating the quality of her experiments:

```
|cs421=> \COPY (SELECT E.NAME, C.DATES, C.FOCUS FROM EXPERIMENT E, CLICASHBOARD C WHERE C.CRID = E.CRID AND C.FOCUS > 3 AND E.NAME = 'Cindy Shaheen' GROUP BY E.NAME, C.FOCUS,C.DATES)TO results.csv WITH CSV
```



## 5 Creative Points

### 5.1 Sophisticated GUI

We created a graphical user interface for our database application. Refer to Section 3.

### 5.2 No Ghost Member

We created a trigger to update the end date of a research assistant when his status goes from active to inactive. This will help professors when a student leaves because they only need to change his status as a 'labmember' and an automatic end date will be generated on the 'supervise' table.

```
cs421=> select * from labmember where name='Zhi Zhang';
      name | contactinformation | status
-----+-----+-----+
 Zhi Zhang | (514)-557-3232 | active
(1 row)

cs421=> select * from supervise where rname='Zhi Zhang';
 startdate | enddate |      rname      |    rname
-----+-----+-----+-----+
 2018-05-01 |           | Cindy Shaheen | Zhi Zhang
(1 row)
```

```

cs421=> CREATE OR REPLACE FUNCTION f_leavingMember()
cs421->   RETURNS trigger AS
cs421-> $BODY$
cs421$> BEGIN
cs421$>   IF NEW.status='inactive' AND OLD.status='active' THEN
cs421$>     UPDATE supervise SET enddate=NOW() WHERE raname=OLD.name;
cs421$>   END IF;
cs421$>   RETURN NEW;
cs421$> END;
cs421$> $BODY$ language plpgsql;
CREATE FUNCTION
cs421=>
cs421=> CREATE TRIGGER leavingMember
cs421->   BEFORE UPDATE
cs421->   ON labmember
cs421->   FOR EACH ROW
cs421->   EXECUTE PROCEDURE f_leavingMember();
CREATE TRIGGER
cs421=> update labmember set status='inactive' where name='Zhi Zhang';
UPDATE 1
cs421=> select * from labmember where name='Zhi Zhang';
      name    | contactinformation | status
-----+-----+-----+
  Zhi Zhang | (514)-557-3232      | inactive
(1 row)

cs421=> select * from supervise where raname='Zhi Zhang';
      startdate | enddate | rname | raname
-----+-----+-----+
  2018-05-01 | 2019-03-22 | Cindy Shaheen | Zhi Zhang
(1 row)

```

As we can see, the created trigger automatically changed the enddate of the research assistant when his status as labmember is changed to inactive from active.